




Maleimide

Product Number: 11306-5, 11306-25, 11306-custom

Product Information Bulletin



Product Background

Alternative Names: [Maleimide, 2,5-Pyrroledione](#)

CAS#: [541-59-3](#)

2,5-Pyrroledione (Maleimide) is a compound used in synthetic applications and in protein chemistry. Maleimide can be derivatized in to larger compounds with specific functionality. The carbon-carbon double bond serves as a useful Michael acceptor and is useful in Diels-Alder applications. Maleimide and its derivatives is also known to be a useful reagent for targeting sulfhydryl groups (-SH) in proteins or other biomolecules.

Product Description

Package Size and Product Numbers

5 grams.....11306-5

25 grams.....11306-25

Product Name:[Maleimide](#)

CAS Number:.....[541-59-3](#)

Product Specifications

- Appearance: White Crystalline Solid
- Gas Chromatography: $\geq 99\%$

CovaChem

6260 East Riverside Blvd

Suite 119

Loves Park, IL 61111

Ph. (815) 315-1271

Fax: (815) 315.1272

Email: info@CovaChem.com