


## TCEP-HCl

Product Number: 11303-1, 11303-5, 11303-25

## Product Information Bulletin



### Product Background

TCEP-HCl ([Tris\(2-carboxyethyl\) phosphine](#)) is a reducing agent commonly used to reduce disulfide bonds in protein molecules. TCEP is soluble in aqueous buffers with at a wide range of pH values. TCEP is a potent reducing agent, that is odorless, and suitable for protein and other biological applications. Unlike other reducing agents, TCEP does not contain any thiol residues (-SH) that may interfere with sensitive assays.

Country of Origin: USA

### Product Description

Package Size and Product Numbers

1 gram.....	11303-1
5 grams.....	11303-5
25 grams.....	11303-25

Product Name: Sulfo-NHS

CAS Number:.....[51805-45-9](#)

### Product Specifications

- Appearance: White Solid
- $H^1$ -NMR: Conform to Structure
- Assay  $\geq$  99 %

CovaChem

6260 East Riverside Blvd

Suite 119

Loves Park, IL 61111

Ph. +1(815) 315-1271

Fax: +1(815) 315.1272

Email: [info@CovaChem.com](mailto:info@CovaChem.com)